

## Product Data Sheet

<b>Product:</b>	CEFT Pool
<b>Product code:</b>	PM-CEFT-3
<b>Source:</b>	Clostridium tetani Epstein-Barr virus (HHV-4) Human cytomegalovirus (HHV-5) Influenza A Selected proteins
<b>No. of peptides:</b>	27 peptides
<b>Amount:</b>	15 nmol (approx. 25ug) / peptide for stimulation of $2,5 \times 10^8$ cells
<b>Length:</b>	diverse
<b>Sequence:</b>	see additional product documentation
<b>Purity:</b>	>90% (HPLC-MS)
<b>Appearance:</b>	white to off-white powder
<b>Remark:</b>	The peptides of this product are supplied as freeze-dried trifluoroacetate salts. For research use only.
<b>Storage conditions:</b>	Store at - 20°C
<b>Application protocol:</b>	<a href="https://jpt.com/support-contact/resources/application-protocols/">https://jpt.com/support-contact/resources/application-protocols/</a>

**Use recommended within 6 months**