



Innovative Peptide Solutions

Product Data Sheet

Product PepMix™ Collection EBV (25ug/peptide)

Product Source PM-C-EBV-2

No. of PepMixes 14

Amount Approx. 25 µg / peptide

Source

vial	Protein - ID	PepMix Name	No. of Peptides
1	P03228	PepMix EBV (BARF1)	53
2	Q04360	PepMix EBV (BMLF1)	117
3	P03191	PepMix EBV (BMRF1)	99
4	P03209	PepMix EBV (BRLF1)	149
5	P03206	PepMix EBV (BZLF1)	59
6	Q8AZK7	PepMix EBV (EBNA-LP)	124
7	P03211	PepMix EBV (EBNA1)	158
8	P12978	PepMix EBV (EBNA2)	119
9	P12977	PepMix EBV (EBNA3a)	234
10	Q1HVG4	PepMix EBV (EBNA3b)	279
11	Q69140	PepMix EBV (EBNA3c)	265
12	P03200	PepMix EBV (GP350/GP340)	224
13	P03230	PepMix EBV (LMP1)	94
14	P13285	PepMix EBV (LMP2)	122

Purity (HPLC) >70% (determined at 220 nm; individual peptides)

Remark The peptides of this product are supplied as trifluoroacetate salts

Storage conditions Store at - 20°C

Solubility recommendations Dissolve in a minimum amount of pure DMSO (approx. 40 µl) and dilute with PBS buffer to the final concentration. Please note that the final concentration of DMSO must be below 1 % (v/v) to avoid toxicity in the biological system.

Use recommended within 6 months.