



Innovative Peptide Solutions

Product Data Sheet

Product PepMix™ Collection EBV (2x1 ug/peptide)

Product Code PM-C-EBV-1

No. of pepmixes 14

Amount Approx. 1 µg / peptide

Source	Row	Column tube 1	Column tube 2	Protein ID	PepMix name	No. of peptides
	A	1	7	P03228	PepMix EBV (BARF1)	53
	B	1	7	Q04360	PepMix EBV (BMLF1)	117
	C	1	7	P03191	PepMix EBV (BMRF1)	99
	D	1	7	P03209	PepMix EBV (BRLF1)	149
	E	1	7	P03206	PepMix EBV (BZLF1)	59
	F	1	7	Q8AZK7	PepMix EBV (EBNA-LP)	124
	G	1	7	P03211	PepMix EBV (EBNA1)	158
	H	1	7	P12978	PepMix EBV (EBNA2)	119
	A	2	8	P12977	PepMix EBV (EBNA3a)	234
	B	2	8	Q1HVG4	PepMix EBV (EBNA3b)	279
	C	2	8	Q69140	PepMix EBV (EBNA3c)	265
	D	2	8	P03200	PepMix EBV (GP350/GP340)	224
	E	2	8	P03230	PepMix EBV (LMP1)	94
	F	2	8	P13285	PepMix EBV (LMP2)	122

Purity (HPLC) >70% (determined at 220 nm; individual peptides)

Remark The peptides of this product are supplied as trifluoroacetate salts

Storage conditions Store at - 20°C

Solubility recommendation Dissolve in a minimum amount of pure DMSO (approx. 2 µl) and dilute with PBS buffer to the final concentration. Please note that the final concentration of DMSO must be below 1 % (v/v) to avoid toxicity in the biological system.

Use recommended within 6 months.