



Innovative Peptide Solutions

Product Data Sheet

Product PepMix Collection HPV E7

Product Source PM-C-HPV-E7

No. of PepMixes 7

Amount Approx. 25 µg / peptide

Source

vial	Protein - ID	PepMix Name	No. of Peptides
1	Q9QLP4	PepMix HPV 06 (E7)	22
2	P04020	PepMix HPV 11 (E7)	22
3	P03129	PepMix HPV 16 (E7)	22
4	P06788	PepMix HPV 18 (E7)	24
5	P27290	PepMix HPV 35 (E7)	22
6	P21736	PepMix HPV 45 (E7)	24
7	P26557	PepMix HPV 58 (E7)	22

Purity (HPLC) >70% (determined at 220 nm; individual peptides)

Remark The peptides of this product are supplied as trifluoroacetate salts

Storage conditions Store at - 20°C

Solubility recommendations Dissolve in a minimum amount of pure DMSO (approx. 40 µl) and dilute with PBS buffer to the final concentration. Please note that the final concentration of DMSO must be below 1 % (v/v) to avoid toxicity in the biological system.

Use recommended within 6 months.