

## Product Data Sheet

<b>Product:</b>	PepMix Collection BKV (25ug/peptide)
<b>Product code:</b>	PM-C-BKV-2
<b>Source:</b>	BK Polyomavirus (BKV) Selected proteins - AGP, LTA, STA, VP1, VP2, VP3, VP4
<b>Protein ID:</b>	P03085, P14999, P03082, P14996, P03094, P03094-2, P03094-3
<b>No. of peptides:</b>	482 peptides (15mers with 11aa overlap)
<b>Amount:</b>	25 µg (approx. 15 nmol) / peptide (for stimulation of up to $2,5 \times 10^8$ cells)
<b>Length:</b>	diverse
<b>Sequence:</b>	Collection of 7 vials a 25ug/peptide; 1 vial of each product.  Product Code: PM-C-BKV-202
<b>Purity:</b>	Development Grade: each peptide purified to > 70% (HPLC/MS)
<b>Appearance:</b>	white to off-white lyophilisate
<b>Remark:</b>	The peptides of this product are supplied as freeze-dried trifluoroacetate salts. For research use only.
<b>Storage conditions:</b>	Store at - 20°C.
<b>Application protocol:</b>	<a href="https://jpt.com/support-contact/resources/application-protocols/">https://jpt.com/support-contact/resources/application-protocols/</a>

**Use recommended within 6 months**