

## Product Data Sheet

|                              |   |
|------------------------------|---|
| <b>Product:</b>              | PepMix Collection Papillomavirus  |
| <b>Product code:</b>         | PM-C-HPV-01   |
| <b>Source:</b>               | Human Papillomavirus (HPV)<br>Selected proteins   |
| <b>Protein ID:</b>           | see datasheet   |
| <b>No. of peptides:</b>      | 396 peptides  |
| <b>Amount:</b>               | 25 µg (approx. 15 nmol) / peptide (for stimulation of up to 2,5 x 10 <sup>8</sup> cells)  |
| <b>Length:</b>               | diverse   |
| <b>Sequence:</b>             | 5 different HPV PepMixes + 2 Controls (15mers, overlap 11)  |
| <b>Purity:</b>               | Development Grade: each peptide purified to > 70% (HPLC/MS)   |
| <b>Appearance:</b>           | white to off-white lyophilisate   |
| <b>Remark:</b>               | The peptides of this product are supplied as freeze-dried trifluoroacetate salts. For research use only.  |
| <b>Storage conditions:</b>   | Store at - 20°C   |
| <b>Application protocol:</b> | <a href="https://jpt.com/support-contact/resources/application-protocols/">https://jpt.com/support-contact/resources/application-protocols/</a> |

**Use recommended within 6 months**

