

## Product Data Sheet

<b>Product:</b>	PepMix Collection Prostate Cancer
<b>Product code:</b>	PM-C-PC
<b>Source:</b>	Human Selected proteins
<b>Protein ID:</b>	see datasheet
<b>No. of peptides:</b>	793 peptides
<b>Amount:</b>	25 µg (approx. 15 nmol) / peptide (for stimulation of up to 2,5 x 10 <sup>8</sup> cells)
<b>Length:</b>	diverse
<b>Sequence:</b>	7 different TAA PepMixes + 2 Controls (15mers, overlap 11); Human (CEA), Human (Metalloreductase STEAP1), Human (Myc protein), Human (PAP), Human (Prostate-specific antigen PSA), Human (TCRgamma alternate reading frame protein TARP), Human (WT1/WT33), CEF Pool (extended), Human (Actin)
<b>Purity:</b>	Development Grade: each peptide purified to > 70% (HPLC/MS)
<b>Appearance:</b>	white to off-white lyophilisate
<b>Remark:</b>	The peptides of this product are supplied as freeze-dried trifluoroacetate salts. For research use only.
<b>Storage conditions:</b>	Store at - 20°C.
<b>Application protocol:</b>	<a href="https://jpt.com/support-contact/resources/application-protocols/">https://jpt.com/support-contact/resources/application-protocols/</a>

**Use recommended within 6 months**