

## Product Data Sheet

<b>Product:</b>	PepMix Collection SARS-CoV-2 Wuhan Structural Proteins (>70%)
<b>Product code:</b>	PM-C-SARS2-WSP-2
<b>Source:</b>	SARS-CoV-2 (Severe Acute Respiratory Syndrome-related coronavirus 2) Wuhan strain Selected proteins - structural proteins
<b>No. of peptides:</b>	643 peptides
<b>Amount:</b>	25 ug (approx. 15 nmol) / peptide for stimulation of $2,5 \times 10^8$ cells
<b>Length:</b>	diverse
<b>Sequence:</b>	8 different SARS-CoV-2 PepMixes: SARS-CoV-2 (Spike Glycoprotein) SARS-CoV-2 (S-RBD) SARS-CoV-2 (NCP) SARS-CoV-2 (VEMP) SARS-CoV-2 (VME1) SARS-CoV-2 (AP3A) SARS-CoV-2 (ORF9B) SARS-CoV-2 (Y14)
<b>Purity:</b>	>70% (HPLC-MS)
<b>Appearance:</b>	white to off-white powder
<b>Remark:</b>	The peptides of this product are supplied as freeze-dried trifluoroacetate salts. For research use only.
<b>Storage conditions:</b>	Store at - 20°C.
<b>Application protocol:</b>	<a href="https://jpt.com/support-contact/resources/application-protocols/">https://jpt.com/support-contact/resources/application-protocols/</a>

**Use recommended within 6 months**