

Product Data Sheet

| | |
|------------------------------|--|
| Product: | PepMix Chlamydia trachomatis (MOMP) |
| Label name: | Chl. trachomatis (MOMP) |
| Product code: | PM-CHLTR-MOMP |
| Source: | Chlamydia trachomatis strain D/UW-3/Cx Outer membrane protein - Major outer membrane porin, serovar D |
| Protein ID: | Q46409 |
| Gene ID: | ompA |
| No. of peptides: | 96 peptides (15mers with 11aa overlap) |
| Amount: | 25 µg (approx. 15 nmol) / peptide (for stimulation of up to 2,5 x 10 ⁸ cells) |
| Length: | 393 aa |
| Sequence: | MKKLLKSVLV FAALSSASSL QALPVGNPAAE PSLMIDGILW EGFGGDPCDP CATWCDAISM RVGYYGDFVF DRVLKTDVNK EFQMGAKPTT DTGNSAAPST LTARENPAYG RHMQDAEMFT NAACMALNIW DRFDVFCTLG ATSGYLKGNs ASFNLVGLFG DNENQKTVKA ESVPNMSFDQ SVVELYTDTT FAWSVGARAA LWECGCATLG ASFQYAQSKP KVEELNVLCN AAEFTINKPK GYVGKEFPLD LTAGTDAATG TKDASIDYHE WQASLALSyr LNMFTPYIGV KWSRASFDAD TIRIAQPKSA TAIFDTTTLN PTIAGAGDVK TGAEGQLGDT MQIVSLQLNK MKSRKSCGIA VGTIIVDADK YAVTVETRLI DERAHVNAQ FRF |
| Purity: | Development Grade: each peptide purified to > 70% (HPLC/MS) |
| Appearance: | white to off-white lyophilisate |
| Remark: | The peptides of this product are supplied as freeze-dried trifluoroacetate salts. For research use only. |
| Storage conditions: | Store at - 20°C. |
| Application protocol: | https://jpt.com/support-contact/resources/application-protocols/ |

Use recommended within 6 months

