

## Product Data Sheet

<b>Product:</b>	PepMix™ HPV18 (E2)
<b>Product code:</b>	PM-HPV18-E2
<b>Source:</b>	Human papillomavirus type 18 Regulatory protein
<b>Protein ID:</b>	P06790
<b>Gene ID:</b>	E2
<b>No. of peptides:</b>	89 peptides
<b>Amount:</b>	25 ug (approx. 15 nmol) / peptide (for stimulation of up to 2,5 x 10 <sup>8</sup> cells)
<b>Length:</b>	365 aa
<b>Sequence:</b>	MQTPKETLSE RLSCVQDKII DHYENDSKDI DSQIQYWQLI RWENAIFFAA REHGIQTLNH QVVPAYNISK SKAHKAIELQ MALQGLAQSA YKTEDWTLQD TCEELWNTEP THCFKKGQQT VQVYFDGNKD NCMTYVAWDS VYMTDAGTW DKTATCVSHR GLYYVKEGYN TFYIEFKSEC EKYGNTGTWE VHFGNNVIDC NDSMCSTSDD TVSATQLVKQ LQHTPSPYSS TVSVGTAKTY GQ TSAATRPG HCGLAEKQHC GPNVPLLGAA TPTGNNKRRK LCSGNTTPII HLGDRNSLK CLRYRLRKHS DHYRDISSTW HWTGAGNEKT GILTVTYHSE TQRTKFLNTV AIPDSVQILV GYMTM
<b>Purity:</b>	>70% (HPLC-MS)
<b>Appearance:</b>	white to off-white powder
<b>Remark:</b>	The peptides of this product are supplied as freeze-dried trifluoroacetate salts. For research use only.
<b>Storage conditions:</b>	Store at - 20°C
<b>Application protocol:</b>	<a href="https://jpt.com/support-contact/resources/application-protocols/">https://jpt.com/support-contact/resources/application-protocols/</a>

**Use recommended within 6 months**