

Product Data Sheet

Product:	PepMaster (Positional Scanning Peptide Library)
Product code:	PM-PSPL
Protein ID:	Not applicable
No. of peptides:	181 different mixtures of peptides
Amount:	440 µg peptide per well
Remark:	This pool consists of 181 peptide pools (9-mers) based on the positional scanning combinatorial library (PSCL) concept for T cell epitope identification.
Purity:	crude
Storage conditions:	Store at - 20°C or lower temperature
Solubility	Dissolve in a minimum amount of pure DMSO (approx. 40 µl) and dilute with PBS buffer or medium to
Recommendation:	the final concentration. Please note that the final concentration of DMSO must be below 1% (v/v) to avoid toxicity in the biological system.
Application protocol:	https://jpt.com/support-contact/resources/application-protocols/

Use recommended within 6 months

Not for human use, for research purposes only

