

Product Data Sheet

Product	PepMaster
Product Code	PM-PSPL
Amount	440 µg peptide per well
Purity (HPLC)	crude
No. of peptides	181 different mixtures of peptides
Protein ID	Not applicable
Remark	This pool consists of 181 peptide pools (9-mers) based on the positional scanning combinatorial library (PSCL) concept for T cell epitope identification.
Storage Conditions	Store at - 20°C or lower temperature
Solubility Recommendation	Dissolve in a minimum amount of pure DMSO (approx. 40 µl) and dilute with PBS buffer or medium to the final concentration. Please note that the final concentration of DMSO must be below 1% (v/v) to avoid toxicity in the biological system.

Use recommended within 6 months.

Not for human use, for research purposes only

<http://www.jpt.com>

JPT Peptide Technologies GmbH
Volmerstrasse 5 (UTZ)
12489 Berlin, Germany

T: +49-30-6392-7878
X: +49-30-6392-7888
e-mail: peptide@jpt.com