

Product Data Sheet

| | |
|------------------------------|---|
| Product: | PepMix™ Pan-HHV6 Select |
| Product code: | PM-Pan-HHV6select-1 |
| Source: | Human herpesvirus 6 (HHV6) Selected proteins |
| Protein ID: | n.a. |
| No. of peptides: | 65 peptides |
| Amount: | 25 ug (approx. 15 nmol) / peptide (for stimulation of up to 2,5 x 10 ⁸ cells) |
| Length: | diverse |
| Sequence: | see additional product documentation |
| Purity: | >70% (HPLC-MS) |
| Appearance: | white to off-white powder |
| Remark: | The peptides of this product are supplied as freeze-dried trifluoroacetate salts. For research use only. |
| Storage conditions: | Store at - 20°C |
| Application protocol: | https://jpt.com/support-contact/resources/application-protocols/ |

Use recommended within 6 months