

## Product Data Sheet

<b>Product:</b>	PepMix Pan-Pertussis Select Pool
<b>Label name:</b>	Pan-Pertussis Select
<b>Product code:</b>	PM-Pan-PERTselect
<b>Source:</b>	Bordetella Pertussis Selected proteins
<b>Protein ID:</b>	n.a.
<b>No. of peptides:</b>	125 peptides (Defined HLA class I & II-restricted T-cell epitopes)
<b>Amount:</b>	15 nmol (approx. 25µg) / peptide for stimulation of 2,5 x 10 <sup>8</sup> cells
<b>Length:</b>	diverse
<b>Sequence:</b>	see additional product documentation
<b>Purity:</b>	Development Grade: each peptide purified to > 70% (HPLC/MS)
<b>Appearance:</b>	white to off-white lyophilisate
<b>Remark:</b>	The peptides of this product are supplied as freeze-dried trifluoroacetate salts. For research use only.
<b>Storage conditions:</b>	Store at - 20°C
<b>Application protocol:</b>	<a href="https://jpt.com/support-contact/resources/application-protocols/">https://jpt.com/support-contact/resources/application-protocols/</a>

**Use recommended within 6 months**

