

## Product Data Sheet

<b>Product:</b>	PepMix SARS-CoV-2 (S-RBD P.1 / Gamma)
<b>Product code:</b>	PM-SARS2-RBDMUT03-1
<b>Source:</b>	SARS-CoV-2 (Severe Acute Respiratory Syndrome-Related Coronavirus 2) (Lineage P.1, B.1.1.28.1, Gamma Variant) Receptor Binding Domain of Spike Glycoprotein - contains mutations K0417T, E0484K, N0501Y
<b>Gene ID:</b>	S
<b>No. of peptides:</b>	53 peptides (15mers with 11 aa overlap)
<b>Amount:</b>	25 µg (approx. 15 nmol) / peptide (for stimulation of up to 2,5 x 10 <sup>8</sup> cells)
<b>Length:</b>	223 aa
<b>Sequence:</b>	RVQPTESIVR FPNITNLCPF GEVFNATRFA SVYAWNRRKRI SNCVADYSVL YNSASFSTFK CYGVSPTKLN DLCFTNVYAD SFVIRGDEV R QIAPGQTGTI ADYNYKLPDD FTGCVIAWNS NNLDSKVGGN YNYLYRLFRK SNLKPFERDI STEIYQAGST PCNGVKGFNC YFPLQSYGFQ PTYGVGYQPY RVVVLSPFELL HAPATVCGPK KSTNLVKNKC VNF □
<b>Purity:</b>	Research Plus Grade: each peptide is guaranteed to be main product (HPLC/MS)
<b>Appearance:</b>	white to off-white lyophilisate
<b>Cytotoxicity:</b>	No cytotoxic effects in resazurin-based cell viability assays with Jurkat cells.
<b>Remark:</b>	The peptides of this product are supplied as freeze-dried trifluoroacetate salts. For research use only.
<b>Storage conditions:</b>	Store at - 20°C
<b>Application protocol:</b>	<a href="https://jpt.com/support-contact/resources/application-protocols/">https://jpt.com/support-contact/resources/application-protocols/</a>

Use recommended within 6 months