

## Product Data Sheet

<b>Product:</b>	PepMix™ SARS-CoV-2 (NS8)
<b>Product code:</b>	PM-WCPV-NS8-1
<b>Source:</b>	SARS-CoV-2 (Severe Acute Respiratory Syndrome-related coronavirus 2) Non-structural protein 8
<b>Protein ID:</b>	PODTC8
<b>Gene ID:</b>	NS8
<b>No. of peptides:</b>	28 peptides
<b>Amount:</b>	25 ug (approx. 15 nmol) / peptide (for stimulation of up to 2,5 x 10 <sup>8</sup> cells)
<b>Length:</b>	121 aa
<b>Sequence:</b>	MKFLVFLGII TTVAAFHQEC SLQSTQHQP YVDDPCPIH FYSKWYIRVG ARKSAPLIEL CVDEAGSKSP IQYIDIGNYT VSCLPFTINC QEPKLGSLVV RCSFYEDFLE YHDVRRVLDLDF I
<b>Purity:</b>	crude (major peak by LC-MS is guaranteed to be peptide of interest for each peptide)
<b>Appearance:</b>	white to off-white powder
<b>Remark:</b>	The peptides of this product are supplied as freeze-dried trifluoroacetate salts. For research use only.
<b>Storage conditions:</b>	Store at - 20°C
<b>Application protocol:</b>	<a href="https://jpt.com/support-contact/resources/application-protocols/">https://jpt.com/support-contact/resources/application-protocols/</a>

**Use recommended within 6 months**