

## Product Data Sheet

|                              |   |
|------------------------------|---|
| <b>Product:</b>              | PepMix™ ZIKV (C) Ultra  |
| <b>Product code:</b>         | PM-ZIKV-C   |
| <b>Source:</b>               | ZIKA virus<br>Capsid protein  |
| <b>No. of peptides:</b>      | 48 peptides   |
| <b>Amount:</b>               | 25 ug (approx. 15 nmol) / peptide (for stimulation of up to $2,5 \times 10^8$ cells)  |
| <b>Length:</b>               | diverse   |
| <b>Sequence:</b>             | see additional product documentation  |
| <b>Purity:</b>               | >70% (HPLC-MS)  |
| <b>Appearance:</b>           | white to off-white powder   |
| <b>Remark:</b>               | The peptides of this product are supplied as freeze-dried trifluoroacetate salts. For research use only.  |
| <b>Storage conditions:</b>   | Store at - 20°C   |
| <b>Application protocol:</b> | <a href="https://jpt.com/support-contact/resources/application-protocols/">https://jpt.com/support-contact/resources/application-protocols/</a> |

**Use recommended within 6 months**