

## Product Data Sheet

<b>Product:</b>	PepStar Candida (MP65)
<b>Product code:</b>	RT-CAND-MP65-2
<b>Source:</b>	Candida Albicans (Yeast) Mannoprotein MP65
<b>Protein ID:</b>	Q9HEP1
<b>No. of peptides:</b>	92 peptides (15mers with 11 aa overlap)
<b>Length:</b>	379 aa
<b>Sequence:</b>	MLFKSFVTFV VLANALAAPL AHQHHQHKEE KRAVHVVTIT NVVVVTIGNG DQTTTFAAPS VAADSSVSVS VNTEPPQNHP TTTQDVASAS TYPSSDTGSA ASSSAAAASS SQAGSEPSGG VSGGAKGIT YSPYSDNGGC KSSSQIASEI AQLSGFNVR LYGVDQDQVA AVLIAKTSSQ KIFAGIFDVS SITSGIESLA EAVKSICGSW DDIYTVSIGN ELVNAGSATP SQIKAYVEEG RKALKAAGYT GPVVSVDTFI AVINNPDLCD YSDYMAVNAH AFFDGHVAAE NSGAWVLQOI QRVWTACGGK KNVLITETGW PSRGDSNGVA VPSKSNQQA ISSIKSSCGA SAILFTAFND LWKADGPYNA EKYWGIYSN

**Specifications:** Microarray glass slide, 1"x3" (2.5 x 7.5cm)

**Layout:** Peptides are printed in 21 mini arrays (3 subarrays each). An exemplary view is shown below. For layout details please refer to .gal-file enclosed to shipment or contact customer support.




figure 1: Multiwell array format: 21 mini arrays per slide. Each mini array displays all peptides in triplicates

**Storage conditions:** Store reffridgerated at 4°C

